


CERTIFICATE OF ANALYSIS

Q-Plex™ Porcine
 Product Name: Cytokine - High
 Sensitivity (4-Plex)
 Order No. 21876
 Part No. 119149PC
 Lot No. PCSK170306 18054 001
 Exp Date: 2023-MAY-14



Software Product Code: PCSM170306CY

Reconstitute Calibrator with 500 µL of Sample Diluent

Component	REF Part No.	LOT Lot No.	Exp. Date
Plate	119152PC	PCSP210514	2023-MAY-14
Sample Diluent	119157PC	PCSS220321	2024-MAR-21
Calibrator	119155PC	PCSM170306	2024-JUN-30
Competitor	N/A	N/A	N/A
Detection	119156PC	PCSD210615	2023-JUN-15
Streptavidin HRP	101173GR	QNSN220111	2023-JUL-11
Substrate A	101159GR	QNSA220309	2024-MAR-09
Substrate B+	101120GR	QNSB+220309	2024-MAR-09
Wash Buffer	101158GR	QNSW220422	2024-APR-22

Spot #	Analyte	Units	Cal Conc	ULOQ	LLOQ	LLD
1	IL-1b	pg/mL	1600	N/A	N/A	2.1
2	IL-6	pg/mL	2000	N/A	N/A	2.5
3	IL-8	pg/mL	2300	N/A	N/A	2.9
4	TNFa	pg/mL	1350	N/A	N/A	1.9
5	REF	N/A	N/A	N/A	N/A	N/A

Quality Release: _____

101896GR 2.0

for Research Use Only