


Q-Plex™ ARRAY

CERTIFICATE OF ANALYSIS

Product Name: Q-Plex™ Human
Malaria (5-Plex)
Order No. 21917
Part No. 565949HU
Lot No. HMRK191125 18162 001
Exp Date: 2023-MAR-21
Software HMRM191125-ML2
Product Code:



Reconstitute Calibrator with 350 µL of Calibrator Diluent

Component	Part No.	Lot No.	Exp. Date
Plate	565952HU	HMRP220801	2024-FEB-01
Caibrator Diluent	565928HU	HMRU211115	2023-MAY-15
Calibrator	565955HU	HMRM191125	2024-FEB-18
Competitor	565990HU	HMRC191014	2024-FEB-18
Assay Diluent	565983HU	HMRV211116	2023-MAY-16
Detection	565956HU	HMRD211119	2023-MAY-19
Streptavidin HRP	101173GR	QNSN210921	2023-MAR-21
Substrate A+	105959GR	MOSA+220810	2024-AUG-10
Substrate B	105920GR	MOSB220810	2024-AUG-10
Wash Buffer	101158GR	QNSW220615	2024-JUN-15

Spot #	Analyte	Units	Cal Conc.	ULOQ	LLOQ	LLD
1	HRP2	pg/mL	2800	2800	0.68	0.53
2	LDH-Pan	pg/mL	67000	67000	16.36	7.28
3	LDH-Pv	pg/mL	19200	19200	4.69	4.23
4	LDH-Pf	pg/mL	20800	20800	5.08	2.33
5	*CRP	ng/mL	38000	38000	9.28	7.65
6	REF	N/A	N/A	N/A	N/A	N/A

*Competitive Assay

Quality Release:


for Research Use Only

101896GR 2.0