

101896GR - 4.1



Quansys Biosciences
365 N 600 W Logan, UT 84321
(435) 752-0531


CERTIFICATE OF ANALYSIS

Q-Plex™ Human Custom 3-Plex
Software Product Code: HA1M200518-IFN

REF Part Number: 107749GR
LOT Lot Number: HCUK200518 18550 001
⌚ Expiration Date: 2024-JAN-07

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	101852GR	HCUP221020	N/A	2024-OCT-20
Sample Diluent	101104HU	HQ6S220107	N/A	2024-JAN-07
Calibrator Diluent	N/A	N/A	N/A	N/A
Assay Diluent	N/A	N/A	N/A	N/A
Calibrator 1	102655HU	HA1M200518	650 µL	2024-OCT-24
Calibrator 2	102855HU	HA2M200303	400 µL	2024-OCT-24
Calibrator Additive	N/A	N/A	N/A	N/A
Controls	N/A	N/A	N/A	N/A
Detection	101856GR	HCUD220908	N/A	2024-SEP-08
Streptavidin	101173GR	QNSN220906	N/A	2024-MAR-06
Substrate A	101159GR	QNSA220825	N/A	2024-AUG-25
Substrate B+	101120GR	QNSB+220825	N/A	2024-AUG-25
Wash Buffer	101158GR	QNSW220909	N/A	2024-SEP-09

Spot #	Analyte	Units	Cal Conc	ULOQ	LLOQ	LLD
1	IFN α	pg/mL	380	380	0.52	0.26
2	IFN β	pg/mL	21700	21700	29.77	14.88
3	IFN γ	pg/mL	3990	3990	5.47	2.74
4	REF	N/A	N/A	N/A	N/A	N/A



Q-Plex™ ARRAYS

QC Release: *[Signature]*
QA Release: *[Signature]*

Date: 21 OCT 22
Date: 21 OCT 22

for Research Use Only