

101896GR - 6.0



Quansys Biosciences  
365 N 600 W Logan, UT 84321  
(435) 752-0531


**CERTIFICATE OF ANALYSIS**

Q-Plex™ Human MMP (6-Plex)  
Software Product Code: MMPM220518

**REF** Part Number: 340949HU  
**LOT** Lot Number: MMPK220518 19191 001  
Expiration Date: 2024-JUN-28

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	340952HU	MMPP230411	N/A	2025-APR-11
Sample Diluent	340957HU	MMPS230331	N/A	2025-MAR-31
Calibrator Diluent	N/A	N/A	N/A	N/A
Assay Diluent	340983HU	MMPV230331	N/A	2025-MAR-31
Competitor	N/A	N/A	N/A	N/A
Calibrator	340955HU	MMPM220518	300 µL	2025-APR-05
Calibrator Additive	N/A	N/A	N/A	N/A
Controls	N/A	N/A	N/A	N/A
Detection	340956HU	MMPD230414	mL	2025-APR-14
Streptavidin HRP	101173GR	QNSN221228	N/A	2024-JUN-28
Substrate A	101159GR	QNSA221221	N/A	2024-DEC-21
Substrate B+	101120GR	QNSB+221221	N/A	2024-DEC-21
Wash Buffer	101158GR	QNSW220803	N/A	2024-AUG-03

Spot #	Analyte	Units	Cal Conc	ULOQ	LLOQ	LLD
1	MMP 1	pg/mL	12563	12563	51.70	25.85
2	MMP 2	pg/mL	85197	85197	350.60	175.30
3	MMP 3	pg/mL	4579	4579	18.84	9.42
4	MMP 7	pg/mL	1745	1745	7.18	3.59
5	MMP 9	pg/mL	44914	44914	184.83	92.42
6	MMP 13	pg/mL	32943	32943	135.57	67.78
7	REF	N/A	N/A	N/A	N/A	N/A



Q-Plex™ ARRAYS

QC Release: [Signature]  
QA Release: [Signature]

Date: 20 APR 23  
Date: 20 APR 23

for Research Use Only