

101896GR - 6.0



Quansys Biosciences
365 N 600 W Logan, UT 84321
(435) 752-0531


CERTIFICATE OF ANALYSIS

Q-Plex™ Bovine Cytokine Panel 1 (5-Plex)
Software Product Code: BC1M230308

REF Part Number: 112349BV
LOT Lot Number: BC1K230308 19358 001
Expiration Date: 2024-FEB-28

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	112352BV	BC1P230228	N/A	2024-FEB-28
Sample Diluent	101104HU	HQ6S220909	N/A	2024-SEP-09
Calibrator Diluent	N/A	N/A	N/A	N/A
Assay Diluent	N/A	N/A	N/A	N/A
Competitor	N/A	N/A	N/A	N/A
Calibrator	112355BV	BC1M230308	500 µL	2024-MAY-24
Calibrator Additive	N/A	N/A	N/A	N/A
Controls	N/A	N/A	N/A	N/A
Detection	112356BV	BC1D230307	N/A	2024-MAR-07
Streptavidin HRP	101173GR	QNSN230216	N/A	2024-AUG-16
Substrate A	101159GR	QNSA230302	N/A	2025-MAR-02
Substrate B+	101120GR	QNSB+230302	N/A	2025-MAR-02
Wash Buffer	101158GR	QNSW230314	N/A	2025-MAR-14

Spot #	Analyte	Units	Cal Conc	ULOQ	LLOQ	LLD
1	IFN γ	pg/mL	1604	1604	2.20	1.10
2	IL-4	pg/mL	477	477	0.65	0.33
3	IL-10	pg/mL	5861	1954	8.04	4.02
4	IL-12p40	pg/mL	210911	210911	289.32	144.66
5	TNF α	pg/mL	3121	3121	12.84	6.42
6	REF	N/A	N/A	N/A	N/A	N/A



Q-Plex™ ARRAYS

QC Release: *[Signature]*
QA Release: *[Signature]*
for Research Use Only

Date: 24 MAY 23
Date: 24 MAY 23