


CERTIFICATE OF ANALYSIS

Q-Plex™ Human Enviromental Enteric Dysfunction (11-Plex)
 Software Product Code: HMDM210721

REF Part Number: 569849HU
LOT Lot Number: HMDK210721 19184 001
 Expiration Date: 2024-AUG-16

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	569852HU	HMDP230406	N/A	2024-OCT-06
Sample Diluent	101104HU	HQ6S220909	N/A	2024-SEP-09
Calibrator Diluent	N/A	N/A	N/A	N/A
Assay Diluent	N/A	N/A	N/A	N/A
Competitor	569890HU	HMDC210708	10 mL	2024-NOV-09
Calibrator	569855HU	HMDM210721	300 µL	2024-NOV-09
Calibrator Additive	N/A	N/A	N/A	N/A
Controls	N/A	N/A	N/A	N/A
Detection	569856HU	HMDD230427	N/A	2024-OCT-27
Streptavidin HRP	101173GR	QNSN230216	N/A	2024-AUG-16
Substrate A	101159GR	QNSA230302	N/A	2025-MAR-02
Substrate B2	101128GR	QNSB2230302	N/A	2025-MAR-02
Wash Buffer	101158GR	QNSW230314	N/A	2025-MAR-14

Spot #	Analyte	Units	Cal Conc	ULOQ	LLOQ	LLD
1	Ferritin	µg/L	111	111	0.15	0.080
2	Tg	µg/L	12.40	12.40	0.017	0.0085
3	sTfR	mg/L	75	75	0.10	0.051
4	HRP2	µg/L	0.96	0.96	0.0013	0.00066
5	IGF-1	ng/mL	202	202	0.28	0.14
6	FGF-21	pg/mL	688	688	0.94	0.47
7	IFABP	pg/mL	2145	2145	2.94	1.47
8	REF	N/A	N/A	N/A	N/A	N/A
9	*CRP	mg/L	5.37	5.37	0.022	0.011
10	*CD-14	ng/mL	2170	2170	8.93	4.47
11	*AGP	g/L	0.30	0.30	0.0012	0.00062
12	*RBP4	µmol/L	1.53	1.53	0.0063	0.0031



Q-Plex™ ARRAYS

*Competitive Assay

QC Release: *Sophia Alvarez* Date: 10 MAY 23
 QA Release: *J. W.* Date: 10 MAY 23

for Research Use Only