

101896GR - 6.0



Quansys Biosciences  
365 N 600 W Logan, UT 84321  
(435) 752-0531

### CERTIFICATE OF ANALYSIS

Q-Plex™ Human Custom 10-Plex  
Software Product Code: HA1M200518-H10

**REF** Part Number: 107749GR  
**LOT** Lot Number: HCUK200518 19229 001  
Expiration Date: 2024-APR-22


Reconstitute Calibrator 1 using 750 µL of Sample Diluent

Reconstitute Calibrator 2 using 650 µL of Calibrator 1

Reconstitute Calibrator 3 using 400 µL of Calibrator 2

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	101852GR	HCUP230330	N/A	2025-MAR-30
Sample Diluent	101104HU	HQ6S220909	N/A	2024-SEP-09
Calibrator Diluent	N/A	N/A	N/A	N/A
Assay Diluent	N/A	N/A	N/A	N/A
Calibrator 1	101855GR	HCUM200804	750 µL	2025-APR-17
Calibrator 2	102655HU	HA1M200518	650 µL	2025-APR-17
Calibrator 3	102855HU	HA2M200303	400 µL	2025-APR-17
Controls	N/A	N/A	N/A	N/A
Detection	101856GR	HCUD230404	N/A	2025-APR-04
Streptavidin HRP	101173GR	QNSN221228	N/A	2024-JUN-28
Substrate A	101159GR	QNSA221221	N/A	2024-DEC-21
Substrate B+	101120GR	QNSB+221221	N/A	2024-DEC-21
Wash Buffer	101158GR	QNSW220422	N/A	2024-APR-22

Spot #	Analyte	Units	Cal Conc	ULOQ	LLOQ	LLD
1	ANG-2	pg/mL	8600	8600	11.80	5.90
2	CEA	pg/mL	62600	6956	85.87	42.94
3	Eotaxin-3	pg/mL	3110	3110	4.27	2.13
4	FAS-Ligand	pg/mL	2360	2360	3.24	1.62
5	IL-5	pg/mL	3100	3100	4.25	2.13
6	IL-7	pg/mL	1680	1680	2.30	1.15
7	MCP-1	pg/mL	1900	633	2.61	1.30
8	MIF	pg/mL	297200	297200	407.68	203.84
9	MIP-1α	pg/mL	2400	800	3.29	1.65
10	TIMP-1	pg/mL	14200	14200	58.44	29.22
11	REF	N/A	N/A	N/A	N/A	N/A



Q-Plex™ ARRAYS

QC Release: [Signature]  
QA Release: [Signature]

Date: 07 JUN 23  
Date: 07 JUN 23

for Research Use Only