

101896GR - 7.0



Quansys Biosciences  
365 N 600 W Logan, UT 84321  
(435) 752-0531

### CERTIFICATE OF ANALYSIS

Q-Plex™ Human Diabetes Panel 1 (2-Plex)  
Software Product Code: HD1M230522




**REF** Part Number: 100849HU

**LOT** Lot Number: HD1K230522 19599 001

**EXP** Expiration Date: 2024-MAY-26

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	100852HU	HD1P230511	N/A	2025-MAY-11
Sample Diluent	100957GR	QAQS230517	N/A	2025-MAY-17
Calibrator Diluent	N/A	N/A	N/A	N/A
Assay Diluent	N/A	N/A	N/A	N/A
Competitor	N/A	N/A	N/A	N/A
Calibrator	100855HU	HD1M230522	400 µL	2024-JUL-18
Calibrator Additive	N/A	N/A	N/A	N/A
Controls	N/A	N/A	N/A	N/A
Detection	100856HU	HD1D230526	N/A	2024-MAY-26
Streptavidin HRP	101173GR	QNSN230628	N/A	2024-DEC-28
Substrate A	101159GR	QNSA230602	N/A	2025-JUN-02
Substrate B+	101120GR	QNSB+230602	N/A	2025-JUN-02
Wash Buffer	101158GR	QNSW230615	N/A	2025-JUN-15

Spot #	Analyte	Cal Conc	Units	ULOQ	LLOQ	LLD
1	GAD-65	14.75	U/mL	8.33	0.034	0.017
2	IA-2	21.40	U/mL	7.67	0.095	0.047
3	Anti-HU IgG	N/A	N/A	N/A	N/A	N/A
4	REF	N/A	N/A	N/A	N/A	N/A



Q-Plex™ ARRAYS

QC Release: [Signature]

Date: 21 JUL 23

QA Release: [Signature]

Date: 21 JUL 23

for Research Use Only