

Q-Plex SARS-CoV-2 Human IgG (5-Plex)

Software Product Code: **HCVM200715 19534 001**

S1 Cutoff Value: 7.7
S2 Cutoff Value: 30
NP Cutoff Value: 25


REF	721649HU
LOT	HCVNK200715 19534 001
	2024-JUL-14

Product Information

Component	Part Number	Lot Number	Exp. Date	Recon. Volume
Plate	721652HU	HCVNP230628	2025-JUN-28	N/A
Sample Diluent	691157HU	HCVS230714	2024-JUL-14	N/A
Calibrator	691155HU	HCVM200715	2025-AUG-02	300 µL
High Control	691131HU	HCVM200715-H	2025-AUG-02	300 µL
Low Control	691132HU	HCVM200715-L	2025-AUG-02	300 µL
Detection	691156HU	HCVD230720	2026-JUL-20	N/A
Streptavidin HRP	691173GR	QNSN230628	2024-DEC-28	N/A
Substrate A	691159GR	QNSA230602	2025-JUN-02	N/A
Substrate B+	691120GR	QNSB+230602	2025-JUN-02	N/A
Wash Buffer	691158GR	QNSW230718	2025-JUL-18	N/A

Assay Characteristics

Analyte	Calibrator Conc. (U/mL)	ULOQ	LLOQ	LLD	CV Limit
S1	1024	700	5.50	5.50	20%
S2	910	700	5.50	5.50	20%
NP	1000	100	5.50	5.50	20%



Lyophilized Control Ranges	Analyte	High Control Range (U/mL)	Low Control Range (U/mL)
	S1	56.96 - 116.59	3.22 - 7.49
	S2	160.30 - 315.90	30.54 - 51.31
	NP	39.68 - 92.60	16.96 - 39.58

Q-Plex™ ARRAY

Quality Control:

[Signature]

Date: 13 NOV 23

Quality Assurance:

[Signature]

Date: 13 NOV 23