

101896GR - 7.0



Quansys Biosciences
365 N 600 W Logan, UT 84321
(435) 752-0531

CERTIFICATE OF ANALYSIS

Q-Plex™ Human Custom 15-Plex
Software Product Code: HA1M220712-L15




REF Part Number: 107749GR
LOT Lot Number: HCUK220712 19877 001
Expiration Date: 2024-OCT-13

Reconstitute Calibrator 1 using 600 µL of Sample Diluent
Reconstitute Calibrator 2 using 400 µL of Calibrator 1

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	101852GR	HCUP230918	N/A	2025-SEP-18
Sample Diluent	101104HU	HQ6S230525	N/A	2025-MAY-25
Calibrator Diluent	N/A	N/A	N/A	N/A
Assay Diluent	N/A	N/A	N/A	N/A
Calibrator 1	102655HU	HA1M220712	600 µL	2025-SEP-28
Calibrator 2	102855HU	HA2M200303	400 µL	2025-SEP-28
Calibrator Additive	N/A	N/A	N/A	N/A
Controls	N/A	N/A	N/A	N/A
Detection	101856GR	HCUD230926	N/A	2025-SEP-26
Streptavidin HRP	101173GR	QNSN230628	N/A	2024-DEC-28
Substrate A	101159GR	QNSA221013	N/A	2024-OCT-13
Substrate B+	101120GR	QNSB+221013	N/A	2024-OCT-13
Wash Buffer	101158GR	QNSW230801	N/A	2025-AUG-01

Spot #	Analyte	Cal Conc	Units	ULOQ	LLOQ	LLD
1	IFN γ	9775	pg/mL	9775	13.41	6.70
2	IL-1 α	3940	pg/mL	3940	5.40	2.70
3	IL-1 β	11800	pg/mL	11800	16.19	8.09
4	IL-6	2850	pg/mL	2850	3.91	1.95
5	IL-8	2080	pg/mL	2080	2.85	1.43
6	IL-10	8400	pg/mL	8400	11.52	5.76
7	IL-12p70	2890	pg/mL	2890	3.96	1.98
8	IL-21	235	pg/mL	235	0.97	0.48
9	IP-10	3240	pg/mL	1080	4.44	2.22
10	MCP-1	1720	pg/mL	573	2.36	1.18
11	MCP-2	585	pg/mL	585	7.22	3.61
12	MIP-1 α	2400	pg/mL	800	3.29	1.65
13	MIP-1 β	1140	pg/mL	1140	1.56	0.78
14	RANTES	6310	pg/mL	6310	8.66	4.33
15	TNF α	6020	pg/mL	6020	8.26	4.13
16	REF	N/A	N/A	N/A	N/A	N/A



Q-Plex™ ARRAYS

QC Release: [Signature]

Date: 28 SEP 23

QA Release: [Signature]

Date: 28 SEP 23

for Research Use Only