

101896GR - 7.0



Quansys Biosciences
365 N 600 W Logan, UT 84321
(435) 752-0531

CERTIFICATE OF ANALYSIS


Q-Plex™ Mouse Custom Stripwell 11-Plex
Software Product Code: MHSM220208-CYR11A



REF Part Number: 107949GR
LOT Lot Number: MCSTK220208 19889 001
Expiration Date: 2024-FEB-12

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	112153MS-15P	MCSTP230707	N/A	2024-JUL-07
Sample Diluent	112457MS	MHSS230512	N/A	2024-JUN-30
Calibrator Diluent	N/A	N/A	N/A	N/A
Assay Diluent	N/A	N/A	N/A	N/A
Competitor	N/A	N/A	N/A	N/A
Calibrator	112455MS	MHSM220208	650 µL	2024-SEP-21
Calibrator Additive	N/A	N/A	N/A	N/A
Controls	N/A	N/A	N/A	N/A
Detection	112456MS	MHSD230212	6 mL	2024-FEB-12
Streptavidin HRP	101173GR	QNSN230628	N/A	2024-DEC-28
Substrate A	101159GR	QNSA221013	N/A	2024-OCT-13
Substrate B+	101120GR	QNSB+221013	N/A	2024-OCT-13
Wash Buffer	101158GR	QNSW230801	N/A	2025-AUG-01

Spot #	Analyte	Cal Conc	Units	ULOQ	LLOQ	LLD
1	GM-CSF	360	pg/mL	470	0.64	0.32
2	IFN γ	660	pg/mL	400	0.55	0.27
3	IL-1 β	1890	pg/mL	1980	2.72	1.36
4	IL-2	170	pg/mL	270	0.37	0.19
5	IL-4	150	pg/mL	150	0.21	0.10
6	IL-6	1180	pg/mL	1180	1.62	0.81
7	IL-10	500	pg/mL	500	0.69	0.34
8	IL-12p70	480	pg/mL	480	0.66	0.33
9	IL-17	620	pg/mL	620	0.85	0.43
10	MIP-1 α	590	pg/mL	590	0.81	0.40
11	TNF α	400	pg/mL	400	0.75	0.38
12	REF	N/A	N/A	N/A	N/A	N/A



Q-Plex™ ARRAYS

QC Release: *[Signature]*

Date: 21SEP23

QA Release: *[Signature]*
for Research Use Only

Date: 21SEP23