

101896GR - 7.0



Quansys Biosciences
365 N 600 W Logan, UT 84321
(435) 752-0531

CERTIFICATE OF ANALYSIS

Q-Plex™ Human Cytokine Inflammation Panel 2 (10-Plex)
Software Product Code: HA1M200518-CIP2




REF Part Number: 110733HU
LOT Lot Number: HCI2K200518 20156 001
Expiration Date: 2024-DEC-19

Reconstitute Calibrator 1 using 650 µL of Sample Diluent
Reconstitute Calibrator 2 using 400 µL of Calibrator 1

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	110752HU	HCI2P221219	N/A	2024-DEC-19
Sample Diluent	101104HU	HQ6S230525	N/A	2025-MAY-25
Calibrator Diluent	N/A	N/A	N/A	N/A
Assay Diluent	N/A	N/A	N/A	N/A
Calibrator 1	102655HU	HA1M200518	650 µL	2025-NOV-30
Calibrator 2	102855GR	HA2M200303	400 µL	2025-NOV-30
Calibrator Additive	N/A	N/A	N/A	N/A
Controls	N/A	N/A	N/A	N/A
Detection	101856GR	HCUD221226	6 mL	2024-DEC-26
Streptavidin HRP	101173GR	QNSN230906	N/A	2025-MAR-06
Substrate A	101159GR	QNSA230302	N/A	2025-MAR-02
Substrate B+	101120GR	QNSB+230302	N/A	2025-MAR-02
Wash Buffer	101158GR	QNSW230919	N/A	2025-SEP-19

Spot #	Analyte	Cal Conc	Units	ULOQ	LLOQ	LLD
1	IFNα	380	pg/mL	380	0.52	0.26
2	IFNβ	21700	pg/mL	21700	29.77	14.88
3	IFNγ	8000	pg/mL	8000	10.97	5.49
4	IL-2	1500	pg/mL	1500	2.06	1.03
5	IL-6	3100	pg/mL	3100	4.25	2.13
6	IL-12p70	2800	pg/mL	2800	3.84	1.92
7	IL-15	8700	pg/mL	8700	11.93	5.97
8	IL-17	19100	pg/mL	19100	26.20	13.10
9	TNFα	5400	pg/mL	5400	7.41	3.70
10	TNFβ	3700	pg/mL	3700	5.08	2.54
11	REF	N/A	N/A	N/A	N/A	N/A



Q-Plex™ ARRAYS

QC Release: [Signature]
QA Release: [Signature]

Date: 04DEC23
Date: 04DEC23

for Research Use Only