

101896GR - 7.0



Quansys Biosciences  
365 N 600 W Logan, UT 84321  
(435) 752-0531

**CERTIFICATE OF ANALYSIS**


Q-Plex™ Human Cytokine Inflammation Panel 2 (10-Plex)  
Software Product Code: HA1M220712-IF2



**REF** Part Number: 110733HU  
**Lot** Lot Number: HCI2K220712 20348 001  
Expiration Date: 2024-AUG-16

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	110752HU	HCI2P221219	N/A	2024-DEC-19
Sample Diluent	101104HU	HQ6S230525	N/A	2025-MAY-25
Calibrator Diluent	N/A	N/A	N/A	N/A
Assay Diluent	N/A	N/A	N/A	N/A
Calibrator 1	102655HU	HA1M220712	600 µL	2026-JAN-16
Calibrator 2	102855HU	HA2M200303	400 µL	2026-JAN-16
Calibrator Additive	N/A	N/A	N/A	N/A
Controls	N/A	N/A	N/A	N/A
Detection	101856GR	HCUD221226	6 mL	2024-DEC-26
Streptavidin HRP	101173GR	QNSN230216	N/A	2024-AUG-16
Substrate A	101159GR	QNSA231006	N/A	2025-OCT-06
Substrate B+	101120GR	QNSB+231006	N/A	2025-OCT-06
Wash Buffer	101158GR	QNSW230718	N/A	2025-JUL-18

Spot #	Analyte	Cal Conc	Units	ULOQ	LLOQ	LLD
1	IFNα	385	pg/mL	385	0.53	0.26
2	IFNβ	21700	pg/mL	21700	29.77	14.88
3	IFNγ	9775	pg/mL	9775	13.41	6.70
4	IL-2	1050	pg/mL	1050	4.32	2.16
5	IL-6	2850	pg/mL	2850	3.91	1.95
6	IL-12p70	2890	pg/mL	2890	3.96	1.98
7	IL-15	9800	pg/mL	9800	13.44	6.72
8	IL-17	19500	pg/mL	19500	26.75	13.37
9	TNFα	6020	pg/mL	6020	8.26	4.13
10	TNFβ	4200	pg/mL	4200	5.76	2.88
11	REF	N/A	N/A	N/A	N/A	N/A



Q-Plex™ ARRAYS

QC Release: [Signature] Date: 18 JAN 24  
 QA Release: [Signature] Date: 1/10/24  
 for Research Use Only