

Q-Plex™ ARRAY

Q-Plex™ SARS-CoV-2 Human IgG (5-Plex)

Software Product Code: HCV200715 20317 001

S1 Cutoff Value: 7.70

S2 Cutoff Value: 30

NP Cutoff Value: 25


REF	721649HU
LOT	HCVNK200715 20317 001
	2025-JAN-11

Product Information

Component	Part Number	Lot Number	Exp. Date	Recon. Volume
Plate	721652HU	HCVNP240110	2026-JAN-10	N/A
Sample Diluent	691657HU	HCVS240111	2025-JAN-11	N/A
Calibrator	691655HU	HCV200715	2026-JAN-15	300 µL
High Control	691631HU	HCV200715-H	2026-JAN-15	300 µL
Low Control	691632HU	HCV200715-L	2026-JAN-15	300 µL
Detection	691656HU	HCV230720	2026-JUL-20	N/A
Streptavidin HRP	691173GR	QNSN231025	2025-APR-25	N/A
Substrate A	691159GR	QNSA231006	2025-MAR-13	N/A
Substrate B+	691120GR	QNSB+231006	2025-MAR-13	N/A
Wash Buffer	691158GR	QNSW240103	2026-JAN-03	N/A

Assay Performance

Analyte	Calibrator Conc. (U/mL)	High Control Range (U/mL)	Low Control Range (U/mL)
S1	1024	41.56 - 117.12	2.41 - 6.80
S2	910	123.62 - 348.38	26.22 - 73.90
NP	1000	33.21 - 93.60	6.01 - 16.93



Quality Control:

[Signature]

Date: 15 JAN 24

Quality Assurance:

[Signature]

Date: 15 JAN 24