

101896GR - 7.0

QUANSYS BIOSCIENCES CERTIFICATE OF ANALYSIS


Quansys Biosciences
365 N 600 W Logan, UT 84321
(435) 752-0531

Q-Plex™ Human Cytokine HS (4-Plex)
Software Product Code: HCHM200803-4P

REF Part Number: 112533HU
LOT Lot Number: HCH4K200803 20582 001
Expiration Date: 2024-DEC-28

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	112552HU	HCH4P231103	N/A	2025-NOV-03
Sample Diluent	101104HU	HQ6S230525	N/A	2025-MAY-25
Calibrator Diluent	N/A	N/A	N/A	N/A
Assay Diluent	N/A	N/A	N/A	N/A
Competitor	N/A	N/A	N/A	N/A
Calibrator	112455HU	HCHM200803	750 µL	2026-MAR-19
Calibrator Additive	N/A	N/A	N/A	N/A
Controls	N/A	N/A	N/A	N/A
Detection	112556HU 01	HCHDQ6230326	6 mL	2025-MAR-26
Streptavidin HRP	101173GR	QNSN230628	N/A	2024-DEC-28
Substrate A	101159GR	QNSA231006	N/A	2025-OCT-06
Substrate B+	101120GR	QNSB+231006	N/A	2025-OCT-06
Wash Buffer	101158GR	QNSW230718	N/A	2025-JUL-18

Spot #	Analyte	Cal Conc	Units	ULOQ	LLOQ	LLD
1	IL-4	18	pg/mL	18	0.025	0.012
2	IL-6	240	pg/mL	240	0.33	0.16
3	IL-10	2320	pg/mL	2320	3.18	1.59
4	IFNγ	93	pg/mL	93	0.13	0.064
5	REF	N/A	N/A	N/A	N/A	N/A



Q-Plex™ ARRAYS

QC Release: [Signature]

Date: 19 MAR 24

QA Release: [Signature]

Date: 19 MAR 24

for Research Use Only