

101896GR - 7.0

QUANSYS
 BIOSCIENCES
CERTIFICATE OF ANALYSIS

Quansys Biosciences
 365 N 600 W Logan, UT 84321
 (435) 752-0531


Q-Plex™ Human Custom 7-Plex
 Software Product Code: HA1M220712-UC7

REF Part Number: 107749GR
LOT Lot Number: HCUK220712 20538 001
 Expiration Date: 2025-APR-25

Reconstitute Calibrator 1 using 600 µL of Sample Diluent
Reconstitute Calibrator 2 using 400 µL of Calibrator 1

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	101852GR	HCUD240221	N/A	2026-FEB-21
Sample Diluent	101104HU	HQ6S240125	N/A	2026-JAN-25
Calibrator Diluent	N/A	N/A	N/A	N/A
Assay Diluent	N/A	N/A	N/A	N/A
Calibrator 1	102655HU	HA1M220712	600 µL	2026-MAR-05
Calibrator 2	102855HU	HA2M200303	400 µL	2026-MAR-05
Calibrator Additive	N/A	N/A	N/A	N/A
Controls	N/A	N/A	N/A	N/A
Detection	101856GR	HCUD240226	N/A	2026-FEB-26
Streptavidin HRP	101173GR	QNSN231025	N/A	2025-APR-25
Substrate A	101159GR	QNSA231205	N/A	2025-DEC-05
Substrate B+	101120GR	QNSB+231205	N/A	2025-DEC-05
Wash Buffer	101158GR	QNSW240103	N/A	2026-JAN-03

Spot #	Analyte	Cal Conc	Units	ULOQ	LLOQ	LLD
1	Eotaxin-1	1740	pg/mL	1740	2.39	1.19
2	IL-1Ra	7930	pg/mL	7930	10.88	5.44
3	IL-6	2850	pg/mL	2850	3.91	1.95
4	IL-8	2080	pg/mL	2080	2.85	1.43
5	IL-10	8400	pg/mL	8400	11.52	5.76
6	IP-10	3240	pg/mL	1080	4.44	2.22
7	MCP-1	1720	pg/mL	573	7.08	3.54
8	REF	N/A	N/A	N/A	N/A	N/A



Q-Plex™ ARRAYS

QC Release: [Signature] Date: 06 MAR 24
 QA Release: [Signature] Date: 6 MAR 24
 for Research Use Only