

101896GR - 7.0

QUANSYS
 BIOSCIENCES
CERTIFICATE OF ANALYSIS


Quansys Biosciences
 365 N 600 W Logan, UT 84321
 (435) 752-0531

Q-Plex™ Human Cytokine Panel 2 (4-Plex)
 Software Product Code: HA1M220712-P24P

REF Part Number: 112033HU
 Lot Number: HCP2K220712 20716 001
 Expiration Date: 2025-JUN-12

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	112052HU	HCP2P240410	N/A	2026-APR-10
Sample Diluent	101104HU	HQ6S240125	N/A	2026-JAN-25
Calibrator Diluent	N/A	N/A	N/A	N/A
Assay Diluent	N/A	N/A	N/A	N/A
Competitor	N/A	N/A	N/A	N/A
Calibrator	102655HU	HA1M220712	600 µL	2026-APR-15
Calibrator Additive	N/A	N/A	N/A	N/A
Controls	N/A	N/A	N/A	N/A
Detection	110956HU 01	HCSDQ6230612	N/A	2025-JUN-12
Streptavidin HRP	101173GR	QNSN240103	N/A	2025-JUL-03
Substrate A	101159GR	QNSA231027	N/A	2025-AUG-23
Substrate B+	101120GR	QNSB+231027	N/A	2025-AUG-23
Wash Buffer	101158GR	QNSW240208	N/A	2026-FEB-08

Spot #	Analyte	Cal Conc	Units	ULOQ	LLOQ	LLD
1	IL-1α	3940	pg/mL	3940	5.40	2.70
2	IL-4	1670	pg/mL	1670	2.29	1.15
3	IL-6	2850	pg/mL	2850	3.91	1.95
4	TNFα	6020	pg/mL	6020	8.26	4.13
5	REF	N/A	N/A	N/A	N/A	N/A



Q-Plex™ ARRAYS

QC Release: [Signature]

Date: 24 APR 24

QA Release: [Signature]

Date: 24 APR 24

for Research Use Only