

101896GR - 7.0



Quansys Biosciences  
365 N 600 W Logan, UT 84321  
(435) 752-0531


**CERTIFICATE OF ANALYSIS**

Q-Plex™ Human Pneumococcal IgG (14-Plex)  
Software Product Code: HPNM231127-W14P

**REF** Part Number: 481849HU  
**LOT** Lot Number: HPNK231127 20328 001  
Expiration Date: 2024-OCT-13

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	481852HU	HPNP231127	N/A	2024-NOV-27
Sample Diluent	481157HU	HPNS240401	N/A	2025-APR-01
Calibrator Diluent	N/A	N/A	N/A	N/A
Assay Diluent	N/A	N/A	N/A	N/A
Competitor	N/A	N/A	N/A	N/A
Calibrator	481155HU	HPNM231127	500 µL	2025-APR-16
Calibrator Additive	N/A	N/A	N/A	N/A
Controls	N/A	N/A	N/A	N/A
Detection	481156HU	HPND231128	N/A	2024-NOV-28
Streptavidin HRP	N/A	N/A	N/A	N/A
Substrate A	101159GR	QNSA231205	N/A	2024-OCT-13
Substrate B+	101120GR	QNSB+231205	N/A	2024-OCT-13
Wash Buffer	101158GR	QNSW240208	N/A	2026-FEB-08

Spot #	Analyte	Cal Conc	Units	ULOQ	LLOQ	LLD
1	1	0.85	µg/mL	0.15	0.0011	0.00054
2	3	0.15	µg/mL	0.10	0.00073	0.00036
3	4	0.33	µg/mL	0.058	0.00042	0.00060
4	5	0.75	µg/mL	0.13	0.00096	0.00048
5	6B	0.91	µg/mL	0.16	0.0012	0.00059
6	7F	0.83	µg/mL	0.15	0.0043	0.0021
7	8	1.42	µg/mL	0.25	0.0071	0.0036
8	9N	0.70	µg/mL	0.12	0.0036	0.0018
9	9V	0.64	µg/mL	0.11	0.0031	0.0016
10	12F	0.22	µg/mL	0.038	0.0011	0.00057
11	14	3.80	µg/mL	0.16	0.0049	0.0024
12	18C	0.73	µg/mL	0.13	0.0037	0.0019
13	19F	1.46	µg/mL	0.26	0.0075	0.0037
14	23F	0.60	µg/mL	0.10	0.0030	0.0015
15	CWPS	N/A	N/A	N/A	N/A	N/A
16	REF	N/A	N/A	N/A	N/A	N/A



**Q-Plex™ ARRAYS**

QC Release:   
QA Release:   
for Research Use Only

Date: 17 APR 24  
Date: 17 APR 24