

101896GR - 7.0



Quansys Biosciences  
365 N 600 W Logan, UT 84321  
(435) 752-0531

### CERTIFICATE OF ANALYSIS


Q-Plex™ Human Cytokine Panel 2 (4-Plex)  
Software Product Code: HA1M220712-P24P



**REF** Part Number: 112033HU  
**LOT** Lot Number: HCP2K220712 20825 001  
Expiration Date: 2025-JUN-12

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	112052HU	HCP2P240410	N/A	2026-APR-10
Sample Diluent	101104HU	HQ6S240125	N/A	2026-JAN-25
Calibrator Diluent	N/A	N/A	N/A	N/A
Assay Diluent	N/A	N/A	N/A	N/A
Competitor	N/A	N/A	N/A	N/A
Calibrator	102655HU	HA1M220712	600 µL	2026-MAY-15
Calibrator Additive	N/A	N/A	N/A	N/A
Controls	N/A	N/A	N/A	N/A
Detection	110956HU 01	HCSDQ6230612	mL	2025-JUN-12
Streptavidin HRP	101173GR	QNSN240103	N/A	2025-JUL-03
Substrate A	101159GR	QNSA240405	N/A	2026-MAR-19
Substrate B+	101120GR	QNSB+240405	N/A	2026-MAR-19
Wash Buffer	101158GR	QNSW231024	N/A	2025-OCT-24

Spot #	Analyte	Cal Conc	Units	ULOQ	LLOQ	LLD
1	IL-1α	3940	pg/mL	3940	5.40	2.70
2	IL-4	1670	pg/mL	1670	2.29	1.15
3	IL-6	2850	pg/mL	2850	3.91	1.95
4	TNFα	6020	pg/mL	6020	8.26	4.13
5	REF	N/A	N/A	N/A	N/A	N/A



Q-Plex™ ARRAYS

QC Release: *John King*

Date: 15 MAY 24

QA Release: *Sally Hansen*

Date: 15 May 24

for Research Use Only