

101896GR - 7.0



Quansys Biosciences  
365 N 600 W Logan, UT 84321  
(435) 752-0531

### CERTIFICATE OF ANALYSIS

Q-Plex™ Human Custom Singleplex  
Software Product Code: FHCM190715-CRTSL 

**REF** Part Number: 107749GR-01


**Lot** Lot Number: HCUK190715 21001 001

**Exp** Expiration Date: 2025-APR-30

#### Reconstitute Calibrator using 550 µL of Complete Sample Diluent

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	101852GR	HCUP240618	N/A	2025-JUN-18
Sample Diluent	101122HU	HQ16240430	N/A	2025-APR-30
Calibrator Diluent	N/A	N/A	N/A	N/A
Assay Diluent	N/A	N/A	N/A	N/A
Competitor	332190HU	FHCC240619	N/A	2025-JUN-19
Calibrator	332155HU	FHCM190715	550 µL	2026-JAN-02
Calibrator Additive	N/A	N/A	N/A	N/A
Controls	N/A	N/A	N/A	N/A
Detection	101856GR	HCUD240621	N/A	2025-JUN-21
Streptavidin HRP	101173GR	QNSN240103	N/A	2025-JUL-03
Substrate A	101159GR	QNSA2310720	N/A	2025-JUL-20
Substrate B+	101120GR	QNSB+230720	N/A	2025-JUL-20
Wash Buffer	101158GR	QNSW231206	N/A	2025-DEC-06

Spot #	Analyte	Cal Conc	Units	ULOQ	LLOQ	LLD
1	Cortisol	714	ng/mL	714	0.98	0.49
2	REF	N/A	N/A	N/A	N/A	N/A



Q-Plex™ ARRAYS

QC Release: 

Date: 03 JUL 24

QA Release: 

Date: 03 JUL 24

for Research Use Only