

101896GR - 7.0

# QUANSYS BIOSCIENCES CERTIFICATE OF ANALYSIS


Quansys Biosciences  
365 N 600 W Logan, UT 84321  
(435) 752-0531

Q-Plex™ Human Cytokine Panel 2 (4-Plex)  
Software Product Code: HA1M220712-P24P

**REF** Part Number: 112033HU  
**LOT** Lot Number: HCP2K220712 21146 001  
Expiration Date: 2025-JUN-12

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	112052HU	HCP2P240410	N/A	2026-APR-10
Sample Diluent	101104HU	HQ6S240125	N/A	2026-JAN-25
Calibrator Diluent	N/A	N/A	N/A	N/A
Assay Diluent	N/A	N/A	N/A	N/A
Competitor	N/A	N/A	N/A	N/A
Calibrator	102655HU	HA1M220712	600 µL	2026-AUG-12
Calibrator Additive	N/A	N/A	N/A	N/A
Controls	N/A	N/A	N/A	N/A
Detection	110956HU_01	HCSDQ6230612	N/A	2025-JUN-12
Streptavidin HRP	101173GR	QNSN240508	N/A	2025-SEP-06
Substrate A	101159GR	QNSA240531	N/A	2026-MAR-19
Substrate B+	101120GR	QNSB+240531	N/A	2026-MAR-19
Wash Buffer	101158GR	QNSW231122	N/A	2025-NOV-22

Spot #	Analyte	Cal Conc	Units	ULOQ	LLOQ	LLD
1	IL-1α	3940	pg/mL	3940	5.40	2.70
2	IL-4	1670	pg/mL	1670	2.29	1.15
3	IL-6	2850	pg/mL	2850	3.91	1.95
4	TNFα	6020	pg/mL	6020	8.26	4.13
5	REF	N/A	N/A	N/A	N/A	N/A



Q-Plex™ ARRAYS

QC Release: *[Signature]*

Date: 12 AUG 24

QA Release: *[Signature]*  
for Research Use Only

Date: 12 AUG 24