



101896GR - 8.0

QUANSYS BIOSCIENCES
CERTIFICATE OF ANALYSIS


Quansys Biosciences
 365 N 600 W
 Logan, UT 84321
 (435) 752 - 0531

Q-Plex™ Mouse Cytokine Panel 1 HS (6-Plex)

REF	Part Number: 111349MS	Software Product Code: MHSM230927-P16P 
LOT	Lot Number: MCP1K230927 21260 001	
	Expiration Date: 2025-FEB-12	

Reconstitute Calibrator using 550 µL of Sample Diluent

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	111352MS	MCP1P240212	N/A	2025-FEB-12
Sample Diluent	112457MS	MHSS240213	N/A	2025-JUN-06
Calibrator	112455MS	MHSM230927	550 µL	2025-SEP-05
Detection	112456MS	MCUD240816	6 mL	2025-AUG-16
Streptavidin HRP	101173GR	QNSN240508	N/A	2025-SEP-06
Substrate A	101159GR	QNSA240531	N/A	2026-MAR-19
Substrate B+	101120GR	QNSB+240531	N/A	2026-MAR-19
Wash Buffer	101158GR	QNSW240702	N/A	2026-JUL-02



Spot #	Analyte	Cal Conc	Units	LLD	LLOQ	ULOQ
1	IFN γ	800	pg/mL	0.76	1.52	595
2	IL-1 α	620	pg/mL	0.70	1.39	546
3	IL-1 β	2390	pg/mL	2.13	4.26	1673
4	IL-6	1110	pg/mL	1.04	2.09	819
5	IL-10	520	pg/mL	0.60	1.19	469
6	TNF α	390	pg/mL	0.28	0.55	217
7	REF	N/A	N/A	N/A	N/A	N/A

Q-Plex™ ARRAYS


QC Release:



Date:

20 SEP 24

QA Release:



Date:

20 SEP 24

for Research Use Only