



101896GR - 8.0

QUANSYS BIOSCIENCES CERTIFICATE OF ANALYSIS


Quansys Biosciences
365 N 600 W
Logan, UT 84321
(435) 752-0531

Q-Plex™ Mouse Cytokine Panel 2 HS (4-Plex)

REF	Part Number: 111149MS	Software Product Code: MHSM230927-P24P 
LOT	Lot Number: MCP4K230927 21261 001	
	Expiration Date: 2025-JUN-06	

Reconstitute Calibrator using 550 µL of Sample Diluent

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	111152MS	MCP4P240906	N/A	2025-SEP-06
Sample Diluent	112457MS	MHSS240213	N/A	2025-JUN-06
Calibrator	112455MS	MHSM230927	550 µL	2025-SEP-13
Detection	101856GR	MCUD240816	6 mL	2025-AUG-16
Streptavidin HRP	101173GR	QNSN240508	N/A	2025-SEP-06
Substrate A	101159GR	QNSA240531	N/A	2026-MAR-19
Substrate B+	101120GR	QNSB+240531	N/A	2026-MAR-19
Wash Buffer	101158GR	QNSW240702	N/A	2026-JUL-02



Spot #	Analyte	Cal Conc	Units	LLD	LLOQ	ULOQ
1	IFN γ	800	pg/mL	0.76	1.52	595
2	IL-1 β	2390	pg/mL	2.13	4.26	1673
3	IL-6	1110	pg/mL	1.04	2.09	819
4	TNF α	390	pg/mL	0.28	0.55	217
5	REF	N/A	N/A	N/A	N/A	N/A

Q-Plex™ ARRAYS

QC Release:

[Signature]

Date: 16 SEP 24

QA Release:

[Signature]
for Research Use Only

Date: 16 SEP 24