



101896GR - 8.0




Quansys Biosciences
365 N 600 W
Logan, UT 84321
(435) 752-0531

Q-Plex™ Human Cytokine Screen (16-Plex)

REF	Part Number: 110933HU	Software Product Code: HA1M220712-CY 
LOT	Lot Number: HCSK220712 21446 001	
	Expiration Date: 2025-APR-25	

Reconstitute Calibrator using 600 µL of Sample Diluent

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	110952HU	HCSP240306	N/A	2026-MAR-06
Sample Diluent	101104HU	HQ6S240125	N/A	2026-JAN-25
Calibrator	102655HU	HA1M220712	600 µL	2026-OCT-15
Detection	110956HU_01	HCSDQ6230725	N/A	2025-JUL-25
Streptavidin HRP	101173GR	QNSN231025	N/A	2025-APR-25
Substrate A	101159GR	QNSA231027	N/A	2025-AUG-23
Substrate B+	101120GR	QNSB+231027	N/A	2025-AUG-23
Wash Buffer	101158GR	QNSW240412	N/A	2026-APR-12



Q-Plex™ ARRAYS

Spot #	Analyte	Cal Conc	Units	LLD	LLOQ	ULOQ
1	IL-1α	3940	pg/mL	2.70	5.40	3940
2	IL-1β	11800	pg/mL	8.09	16.19	11800
3	IL-2	5150	pg/mL	3.53	7.06	5150
4	IL-4	1670	pg/mL	1.15	2.29	1670
5	IL-5	2950	pg/mL	2.02	4.05	2950
6	IL-6	2850	pg/mL	1.95	3.91	2850
7	IL-8	2080	pg/mL	1.43	2.85	2080
8	IL-10	8400	pg/mL	5.76	11.52	8400
9	IL-12p70	2890	pg/mL	1.98	3.96	2890
10	IL-13	2865	pg/mL	1.97	3.93	2865
11	IL-15	9800	pg/mL	6.72	13.44	9800
12	IL-17	19500	pg/mL	13.37	26.75	19500
13	IL-23	45950	pg/mL	31.52	63.03	45950
14	IFNγ	9775	pg/mL	6.70	13.41	9775
15	TNFα	6020	pg/mL	4.13	8.26	6020
16	TNFβ	4200	pg/mL	2.88	5.76	4200
17	REF	N/A	N/A	N/A	N/A	N/A

QC Release: 
 QA Release: 
 for Research Use Only

Date: 16 OCT 24
 Date: 16 OCT 24