



101896GR - 8.0

QUANSYS BIOSCIENCES
CERTIFICATE OF ANALYSIS


Quansys Biosciences
 365 N 600 W
 Logan, UT 84321
 (435) 752-0531

Q-Plex™ Mouse Cytokine Panel 2 HS (4-Plex)

REF	Part Number: 111149MS	Software Product Code: MHSM230927-P24P 
LOT	Lot Number: MCP4K230927 21571 001	
	Expiration Date: 2025-JUN-06	



Reconstitute Calibrator using 550 µL of Sample Diluent

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	111152MS	MCP4P240906	N/A	2025-SEP-06
Sample Diluent	112457MS	MHSS240213	N/A	2025-JUN-06
Calibrator	112455MS	MHSM230927	550 µL	2025-NOV-13
Detection	112456MS	MHSD240919	6 mL	2025-SEP-19
Streptavidin HRP	101173GR	QNSN240103	N/A	2025-JUL-03
Substrate A	101159GR	QNSA231027	N/A	2025-AUG-23
Substrate B+	101120GR	QNSB+231027	N/A	2025-AUG-23
Wash Buffer	101158GR	QNSW231218	N/A	2025-DEC-18



Spot #	Analyte	Cal Conc	Units	LLD	LLOQ	ULOQ
1	IFN γ	800	pg/mL	0.55	1.10	800
2	IL-1 β	2390	pg/mL	1.64	3.28	2390
3	IL-6	1110	pg/mL	0.76	1.52	1110
4	TNF α	390	pg/mL	0.27	0.53	390
5	REF	N/A	N/A	N/A	N/A	N/A

Q-Plex™ ARRAYS

QC Release:  Date: 13 NOV 24
 QA Release:  Date: 13 NOV 24
 for Research Use Only