

101896GR - 8.0




Quansys Biosciences
365 N 600 W
Logan, UT 84321
(435) 752-0531

Q-Plex™ Human Cytokine Panel 1 (6-Plex)

REF	Part Number: 115233HU	Software Product Code: HA1M220712-P16P
LOT	Lot Number: HCP1K220712 21616 001	
	Expiration Date: 2025-JUL-25	

Reconstitute Calibrator using 600 µL of Sample Diluent

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	115252HU	HCP1P231207	N/A	2025-DEC-07
Sample Diluent	101104HU	HQ6S240805	N/A	2026-AUG-05
Calibrator	102655HU	HA1M220712	600 µL	2026-DEC-04
Detection	110956HU_01	HCSDQ6230725	N/A	2025-JUL-25
Streptavidin HRP	101173GR	QNSN241105	N/A	2026-MAY-05
Substrate A	101159GR	QNSA240801	N/A	2026-MAR-19
Substrate B+	101120GR	QNSB+240801	N/A	2026-MAR-19
Wash Buffer	101158GR	QNSW240809	N/A	2026-AUG-09



Spot #	Analyte	Cal Conc	Units	LLD	LLOQ	ULOQ
1	IL-1β	11800	pg/mL	8.09	16.19	11800
2	IL-6	2850	pg/mL	1.95	3.91	2850
3	TNFα	6020	pg/mL	4.13	8.26	6020
4	IFNγ	9775	pg/mL	6.70	13.41	9775
5	IL-1α	3940	pg/mL	2.70	5.40	3940
6	IL-10	8400	pg/mL	5.76	11.52	8400
7	REF	N/A	N/A	N/A	N/A	N/A

Q-Plex™ ARRAYS

QC Release: Date: 04 DEC 24
 QA Release: Date: 04 DEC 24
for Research Use Only